ERNST HINRICH'S Dental GmbH

Borsigstraße 1 · DE-38644 Goslar/Germany
Tel. +49 (0) 53 21 / 5 06 24-25 · Fax +49 (0) 53 21 / 5 08 81
info@hinrichs-dental.de · www.hinrichs-dental.de

Hinrvest® RP

DIN EN ISO 15912, Typ 2, Klasse 2

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DIN EN ISO 15912, type 2, class 2

Instruction for use

A super fine phosphate bonded, graphite free precision investment for the production of partial denture frameworks. Hinriwest® RP is developed especially for the casting of printed (rapid prototyping) and milled wax & resin parts.

Hinriwest® RP is usable for the conventional and speed casting techniques and is suitable for casting all types of dental alloys (except Titanium).

Technical data:

Mixing liquid:	Hinriwest® Liquid
Mixing ratio, powder : liquid:	100 g : 20 ml
Mixing under vacuum:	60 sec.
Working time:	4 - 5 min.
Working temperature:	approx. 22 °C (Room temperature)
Initial setting time:	7 - 10 min.
Setting expansion:	approx. 1,00 % (75% diluted type 100 liquid)
Preheating temperature range for placing in the furnace:	900 - 950 °C
Place in the hot furnace:	20 min. after investing

Concentration of the expansion liquid: 20 ml liquid / 100 g powder

Frameworks - Liquid ratio Typ 100 / distilled water: 80 % / 20 %

100 g	180 g	400 g
16 ml / 4 ml	29 ml / 7 ml	64 ml / 16 ml

Investing: Follow the instructions according to the SILADENT-system when attaching sprues to the framework. Place the framework in the model forming silicone sleeves (REF 892003 - 892005) and ensure that it is placed out of the heat centre of the mould. Thin sections, e.g. clasps, should face towards the sides of the mould and be 5 mm from the casting ring. Measure and mix Hinriwest® RP as described above. Pour with minimum vibration until the framework is completely covered. Switch off the vibrator, fill the ring and do not vibrate again. **Caution: A pressure-free investing is important!**

1888

Preheating with the speed casting technique: **20 minutes** after investing, place the mould with the sprue hole pointing downwards in a furnace preheated to the correct temperature of **900° C – 950° C**. ERNST HINRICHSEN Dental recommend a maximum preheating temperature of **950° C**. This temperature is perfectly adequate for both casting the mould and final thermal expansion. Moulds are ready to cast approx. 60 minutes after the final temperature is reached. During the first 15 minutes after inserting the moulds, the furnace fan and extractor should remain switched off and the furnace door should not be opened because of the risk of instant combustion.

Please note: Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. ERNST HINRICHSEN Dental products are subjected to strict quality controls. We reserve the right to make technical changes.

Warning: Investments contain quartz and cristobalite. Do not breathe in dust. Danger of lung diseases (silicosis or cancer). Use a suitable dust mask!

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